





【開始】 最初に読んでください 	[2] 	[1] 
サージスクリーン 3SS851 Exp:2026/08/04 	[3] 	ディエコ A血球 DIA736 Exp:2026/08/04 
[1] 	サージスクリーン 3SS857 Exp:2026/08/04 	[1] 
[2] 	[1] 	ディエコ A血球 DIA737 Exp:2026/08/04 
[3] 	[2] 	[1] 
サージスクリーン 3SS851Z Exp:2026/08/04 	[3] 	リゾルブパネル B RB710 Exp:2026/08/04 
[1] 	バイオビュースクリーン BS4531 Exp:2026/08/03 	[12] 
[2] 	[1] 	[13] 
[3] 	[2] 	[14] 
サージスクリーン 3SS852 Exp:2026/08/04 	[3] 	[15] 
[1] 	バイオビュースクリーン BS4532 Exp:2026/08/03 	[16] 
[2] 	[1] 	[17] 
[3] 	[2] 	[18] 
サージスクリーン 3SS852Z Exp:2026/08/04 	[3] 	[19] 
[1] 	バイオビュースクリーン BS4533 Exp:2026/08/03 	[20] 
[2] 	[1] 	[21] 
[3] 	[2] 	[22] 
サージスクリーン 3SS855 Exp:2026/08/04 	[3] 	セレクトシエン S641 Exp:2026/08/04 
[1] 	バイオビュースクリーン BS4536 Exp:2026/08/04 	[1] 
[2] 	[1] 	[2] 
[3] 	[2] 	【終了】 最後に読んでください 
サージスクリーン 3SS856 Exp:2026/08/04 	[3] 	
[1] 	ディエコ A血球 DIA735 Exp:2026/08/04 	

※このバーコードで登録したLotをQuickAntigamで「血球試薬の追加」をする際は、DonorNoのアルファベットを数字の9に置き換えて入力してください。

※0.8%リゾルブパネルAおよびCは通常のリゾルブパネルAおよびCとして、アブテストセル皿はサージスクリーンとして読み込まれます。